

Infinite Insight XCUBE 90 on GPX



ALPINION MEDICAL SYSTEMS

We are Ultrasound Professionals






X-CUBE 90 on GPX offers the highest level of efficiency and insight for healthcare providers, enabling you to utilize clear and detailed image and a comfortable user interface.

X-CUBE 90 on GPX reflects the demands of fields that require quick and accurate diagnosis. With ALPIN-ION's amassed image processing technology and the latest diagnostic solution, more precise diagnosis and quicker decision are now possible.

Improve the quality of your diagnosis with X-CUBE 90 on GPX. ALPINION



eXtensive Insight

Powered by X⁺ Architecture

X⁺ Architecture is a premium imaging platform created by combining

X⁺ Crystal Signature[™], the transducer technology with high sensitivity, wideband, and


X⁺ FIT, ALPINION's cutting edge beamforming and data processing technology.

The X-CUBE 90 on GPX has undergone a comprehensive system upgrade that integrates advanced image processing technology X⁺ Focus.

X⁺ Focus reduces noise and artifacts to provide excellent contrast resolution delicately.

The most important factor for accurate and confident diagnosis is image quality.

X⁺ Architecture supports your expertise.



X⁺ Crystal Signature transducers generate high sensitive, wide-band signals

X⁺ FIT, high performance beamforming engine to generate high-resolution images from massive scan data



X⁺ Focus, advanced image processing technology provides clearer images in a wide range of different diseases







Resolution


Penetration

Uniformity

X⁺ Crystal Signature™




Upgraded single crystal and a special material matching layer and backing material allow heat dissipation and minimalize signal loss. ALPINION's signature transducer technology allows you to experience insightly clear image quality with high sensitivity and wideband, which shows significant improvements in image resolution, contrast, and uniformity. ALPINION develops and manufactures its own transducers, the key to ultrasound imaging.








 X^{+} FIT is a beamforming technology which improves resolution, frame rate, and uniformity. X^{+} FIT utilizes massive parallel beamforming whose architecture is capable of transmitting 4 times larger volume with 10 times faster capacity and also processing data 14 times faster than our previous model (E-CUBE 15 Platinum). X-CUBE 90 on GPX's high resolution imaging is a prerequisite for accurate diagnoses and quick decision-making.




X⁺ Focus


 X^{+} Focus is new image processing technology that displays a clearer image with massive scan data. It realizes optimal image quality by enhancing true clinical information and suppressing unnecessary noise and other artifacts from rich scan data. In particular, X^{+} Focus contributes to accurate diagnosis by providing adequate information, well-defined border, sharp contrast, and rich gray-scale without compromising on maintained tissue information for diagnosis.



Powered by X⁺ Architecture



Hepatic Vein in CF mode




X⁺ MicroView of Portal Vein




Gallbladder in B mode



Pancrease in B mode




CEUS with Liver



Breast Nodule in B mode




Mitral Valve Short Axis




Tricuspid Valve Regusitation




Carotid Artery




Brilliant Flow of Thyroid Nodule




Biceps Tendon




Quadriceps Tendon





Shear Wave Elastography (SWE)

Shear Wave Elastography evaluates tissue stiffness. This non-invasive and advanced tool provides quantitative assessment for chronic liver disease and suspicious nodule of Thyroid and Breast. The reliable and reproducible result of SWE helps to increase diagnostic confidence.



Attenuation Tissue Imaging (ATI)


The ATI quantifies the attenuation to measure the level of steatosis with color-mapped image. The ATI provides quickly calculated value with high reproducibility for diagnosis efficiency. Nonalcoholic fatty liver(NAFLD) is the most common cause of liver dysfunction, with a prevalence of up to 20% of the total population in world wide.




Contrast-Enhanced Ultrasound (CEUS) with TIC

This is a function to diagnose patients using various angiographic patterns that appear while a contrast medium, administered intravenously, diffuses in blood vessels and organ tissue. CEUS has many advantages in various clinical indications for liver disease.


In addition, Time Intensive Curve(TIC)-analysis available with CEUS is a quantitative tool for the dynamic evaluation of target lesions.






X⁺ MicroView

 X^* MicroView is the vascular imaging mode which displays micro blood flow. Users can observe the low-speed blood flow of tiny blood vessel. With this technology, low-speed blood flow areas which users previously couldn't see on Color Doppler are shown at the fast frame rate.



Brilliant Flow

It is an special visualization impact technology of blood flow. It provides the intuitive understanding of vessel structure and boundaries by images of blood vessels that look like 3D on a 2D gray scale image. It may provide additional spatial information in CF, PDI or X⁺ MicroView.



CUBE Strain[™]

This is a non-invasive examination method that is used to assess the myocardial function more objectively. Users can track speckles in 2D heart images, digitize the movement of each myocardial segment, and check quantified data.



X⁺ Assistant

 X^{+} Assistant enables users to reduce keystrokes by at least 50% and save time when conducting examinations. Registered optimal scan protocols that follow guidelines that each application requires help to reduce operator's fatigue and increase efficiency in operation. In addition, it is also possible to register applications and protocols optimized.

X⁺ Compare

This is a feature that allows users to import patients' previous study from a PACS server or hard disk and compare when scanning and reviewing. By comparing previous images of the same region, users can easily observe the patient's progress, thereby making this feature useful for providing patient care. X* Compare can be used not only in scan mode, but in review mode (E-View) as well. * X* Compare supports ultrasound studies only.

Auto Measurement Solutions

Such as Auto IMTTM, Auto EF, and an assessment system for the breast and thyroid image enhance productivity.

Needle Vision[™] Plus

This feature shows the shape and orientation of the needle clearly in real time. Using Beam Steering technology, detect and locate the needle precisely.

USB Real-time

USB real-time recording makes data storage easier by allowing users to record ultrasound scan images on USB memory in real time. Videos are recorded as high-definition and stored in system quickly.

eXceptional Care

The detailed work process designed by X-CUBE 90 on GPX cares about realizing a more efficient working environment for medical professionals. Comprehensive tools such as auto measurement, analytics, exam assists, and ergonomic design increase productivity across a wide range of patient cases and applications.

Power Preset

With a single touch, users can load a transducer with allocated image preset. Quick and easy application of preset shorten the image setup time.

Time Saving Solution

Quick system ready by fast boot-up in 1 minute. And the built-in battery guarantees mobility while it is in use. When not using the system, Sleep mode helps save energy and quickly boot up in 6sec.

Antivirus Solution

As X-CUBE 90 on GPX's automatic virus scanning function is activated whenever the system boots, patients' information is more securely protected.

On-line Connection

With the Online Connection solution, direct support for customers through online chat and face-to-face talk is possible anywhere.

Wide FOV

The wide field of view provides more anatomical structures in real time. The extended view helps better understanding of the relationship between organs and additional anatomical information improves value of a report.



Intuitive design for your daily practice





23 inch FHD LED Screen

Enjoy a larger widescreen and clear, highdefinition images.



Compact Size & Weight

Its compact size and light weight measuring 85 kg provides great mobility.



Motorized Control Panel

The control panel can be easily adjusted vertically with its motorized button.



0 00

4 Swivel Wheels

Thanks to the 4 swivel wheels, the system can be moved around with a little effort and can be operated comfortably in any position.



Additional Probe Holder

Providing more ease to pick up the transducer and gel, the additional holder provides a convenient diagnostic environment to fit the body shape of user.

12.1 inch Tilting LCD Color Touch Panel

We've combined angle adjustment functionality with its larger, full-color touch panel to maximize user convenience.

320° Articulating Monitor Arm & ±90°, 60° Angle Adjustable Gel Warmer


Ergonomic design helps to promote healthy posture, which reduces musculoskeletal pain.

5 Transducer Ports

You can use a variety of transducers as required with no need to replace them, and thus saving time during conducting diagnosis.

ALPINON'sTransducer Technology

The key to transducer technology is its process of generating a high sensitivity and wideband acoustic signal and receiving it without any loss, then converting it to a digital signal. ALPINION MEDICAL SYSTEMS has researched and developed the core technology of transducers and X⁺ Crystal Signature™ is the accomplishment of all the innovative technologies we've accumulated.



Transducers H61-E7S Developed and manufactured by ALPINION





SC1-7H

X⁺ Crystal Signature™

High Density Single Crystal Convex Abdomen, OB/GYN, Pediatric, EM, Urology



SC2-9H No



High Density Single Crystal Convex Abdomen, OB/GYN, Pediatric, EM, Urology



SC2-11H Ne

X⁺ Crystal Signature™

High Density Single Crystal Microconvex

Abdomen, OB/GYN, Pediatric, EM

Volume Convex



SVC1-8H

Single Crystal Volume Convex Abdomen, OB/GYN, Pediatric, EM, Urology

Linear



SL3-19X New

X⁺ Crystal Signature™ Extreme High Density Single

Crystal Linear, Seamless window Abdomen, Pediatric, OB/GYN, EM, MSK, Vascular, Small Parts,



SL3-19H New

X⁺ Crystal Signature™ (5) Extreme High Density Single Crystal Linear, Seamless window

Abdomen, Pediatric, OB/GYN, EM, MSK, Vascular, Small Parts,



L3-12X

Hgh Density Linear Pediatric, EM, MSK, Vascular, Small Parts



L3-8H

Hgh Density Low Freq. Linear Breast, EM, MSK, Vascular, Small Parts



L10-25H

Wideband Ultra High Crystal Linear, Seamless window

EM. MSK. Vascular. Small Parts

Endocavity



EC2-11H New

X⁺ Crystal Signature™

High Density Single Crystal Endocavity(Straight), Small tip, FOV MAX.230°

OB/GYN, Urology, EM



EV2-11H New

X⁺ Crystal Signature™

High Density Single Crystal Endocavity(Curved), Small tip, FOV MAX.230°

OB/GYN, Urology, EM



VE3-10H New

Volume Endocavity OB/GYN, Urology, EM

Phased Array



MP1-5X New



X⁺ Crystal Signature™

Abdomen, Pediatric, Cardiac, EM, TCD



SP3-8T

Single Crystal Phased Array Abdomen, Cardiac, EM, Pediatric

Pencil



CW2.0

Pencil Typed Cardiac

CW5.0

Pencil Typed Cardiac



CW8.0

Pencil Typed





XCUBE 90 on GPX

eXtensive Insight eXpanded Capability eXceptional Care

ALPINION MEDICAL SYSTEMS Co., Ltd.

15, Magokjungang 14-ro, Ganseo-gu, Seoul, Korea

Homepage www.alpinion.com

E-mail international@alpinion.com

TEL +82-2-3777-8600 FAX +82-2-3777-8691

Standalone clinical images may have been cropped to better visualize pathology.

Copyright@2022 ALPINION MEDICAL SYSTEMS CO., LTD. All rights reserved.

Catalogue contents may change without prior notice to customers due to performance enhancements.

